


UV-Vis Spectrophotometer

Perkin Elmer Lambda 25



Double-beam scanning UV-Vis Spectrophotometer, PC operated via UVWinLab software (manufacturer Perkin Elmer, 2006).

- Wavelength range: 190-1100 nm
- Wavelength precision: ± 0.1 nm
- Lamp change: automatic at 326 nm
- The scanning speed: can vary from 7.5 to 2880 nm/min
- Photometric precision: Absorbance ± 0.003



It is used for the qualitative identification of a compound in a solution, as well as for the quantitative identification of a substance in a solution.

The principle of the UV-Vis spectroscopy is based on the Lambert-Beer Law, which states that there is a linear relationship between the concentration and the absorbance of the solution, which enables the concentration of a solution to be calculated by measuring its absorbance.

LAMBERT-BEER LAW

T= transmission

$$T = \frac{I}{I_0} \times 100$$

A= absorbance

$$A = \frac{I_0 - I}{I_0} \times 100$$

E= extinction (molar absorptivity) $E = \log \frac{I_0}{I} = \epsilon lc$

c = the concentration of the irradiated solution (mol/L)

l = length of light path (cm)

ϵ = molar absorptivity or molar extinction coefficient



Dr. Lilia Anghel
Office: 220
Email: lilia.anghel@ichem.md

Institute of Chemistry,
str. Academiei, 3,
Chisinau, Republic of Moldova
Web: www.ichem.md